



University
of Glasgow

Imperial College
London



Visualizing Dynamic Magnetism in Nanostructures using Electron Microscopy

Trevor P. Almeida

Materials and Condensed Matter Physics Group

School of Physics and Astronomy, University of Glasgow, G12 8QQ, UK

EPSRC





DENS
solutions



**NATURAL
ENVIRONMENT
RESEARCH COUNCIL**


- **Motivation**
- **Transmission Electron Microscopy**
- **Magnetic Minerals**
 - **Electron holography**
- **Phase transition in FeRh thin films**
 - **Differential phase contrast imaging (DPC)**




The demand for **improved functionality** and **reduction in size** of a range modern devices has led to the rapid development of **new magnetic materials**

→ driving the need to **visualise localised magnetism on the nanoscale**


Magnetic data storage




Nanomedicine




Earth science





Example TEM schematic

One of many types of TEMs



Sample thickness
< 200 nm



Resolution
~ 80 pm


Information:

- Morphology
- Chemistry
- Structure
- Magnetism
- Electronic
- Etc.




- Motivation
- Transmission Electron Microscopy
- **Magnetic Minerals**
 - Electron holography
- Phase transition in FeRh thin films
 - Differential phase contrast imaging (DPC)





between reversals




during a reversal

To interpret palaeomagnetic data we need to understand mechanisms that:


1. induce magnetic remanence
2. alter magnetic remanence


1. Chemical alteration –
change in oxidation state

2. Thermomagnetic behaviour



after Butler (1982)






- FEI Titan Analytical
- C_s correction on condenser lens, *i.e.* probe corrected
- Operated at 300kV
- HR-STEM and chemical mapping at atomic scale
- Biprism and Lorentz lens for electron holography of magnetic fields

DTU




DTU Cen
Center for Electron Nanoscopy




Phase shift: $\varphi(X) = C_E \int V(x, z) dz - \left(\frac{e}{h}\right) \iint B_{\perp}(x, z) dx dz$

Mean inner potential
Magnetic induction


Total phase shift



Mean inner potential




Magnetic contribution





—

=

Vortex state and stray magnetic field – PSD particle



Oxidation of Fe_3O_4 :





● A sites
tetrahedral
Fe³⁺

● B sites
octahedral
Fe³⁺, Fe²⁺


Ferrimagnetic




Oxidation of Fe_3O_4 :




Slightly weaker ferrimagnetic



Oxidation of Fe_3O_4 :



Canted antiferromagnetic






- FEI Titan E-Cell
- C_s correction on objective lens, *i.e.* image corrected
- Operated at 300kV
- Various gases, *i.e.* H_2 , He and H_2O up to 1000 Pa
- Heating specimen holder < 1000°C




DTU Cen
Center for Electron Nanoscopy


Energy dispersion of 0.02eV & resolution of 0.5eV



EEL spectrum from
Fe₃O₄ particles




EEL spectrum from
γ-Fe₂O₃ sample




EEL spectrum from
reference **α-Fe₂O₃** sample

Heated at **700°C**, under **9mbar O₂** for **8 hours**





Heated at **700°C**, under **9mbar O₂** for **8 hours**



Before oxidation

Oxidised








- FEI Titan HOLO
- C_s correction on objective lens, *i.e.* image corrected
- Operated at 60 - 300kV
- 3 biprisms and Lorentz lens for electron holography of magnetic fields
- 11 mm pole piece gap to allow tilting to $\pm 75^\circ$
- Specifically designed to allow for *in situ* experiments

Ernst Ruska-Centrum
für Mikroskopie
und Spektroskopie
mit Elektronen

ER-C




Wildfire heating holder and EMheaterchips™.




Ultra stable holder during heating up to $< 1000^{\circ}\text{C}$

TRM = total phase shift – MIP (measured at each temperature in separate experiment)


Fe_3O_4 grain at remanence upon heating to **550°C** and cooling to room temperature





Fe_3O_4 grain at remanence upon heating to **600°C** and cooling to room temperature





Combining electron holography with liquid-cell TEM holder





Combining electron holography with liquid-cell TEM holder



- Motivation
- Transmission Electron Microscopy
- Magnetic Minerals
 - Electron holography
- **Phase transition in FeRh thin films**
 - Differential phase contrast imaging (DPC)



Equiatomic intermetallic iron-rhodium compound ($\text{Fe}_{48}\text{Rh}_{52}$ to $\text{Fe}_{56}\text{Rh}_{44}$)



380 K: AF \rightarrow FM

Fe: $\sim 3.1 \mu_B$

Rh: $\sim 1 \mu_B$

Zakharov A I et al., *J. Exp. Theor. Phys.*, 1964
 Kouvel J S et al., *J. Appl. Phys.* 1962



DC magnetron sputter co-deposition

- FeRh targets
- MgO and GaAs substrates
- Samples:
 - 1) FeRh on MgO
 - 2) Planar FeRh TEM samples via HF-etching




UNIVERSITY OF LEEDS







Focused Ion Beam / Scanning Electron Microscope (FIB-SEM)

- Cross-sectional and planar FeRh samples



Lightning holder with heater / biasing MEMS chip






- JEOL ARM200cF - "MagTEM"
- C_s correction on condenser lens, *i.e.* probe corrected
- Operated at 60kV - 200kV
- HR-STEM and chemical mapping at atomic scale, with EDX and EELS
- Lorentz lens and segmented / pixelated detectors for imaging of magnetisation



HF-etched FeRh planar sample as a function of temperature

20°C \longrightarrow 140°C \longrightarrow 20°C


Bright field TEM imaging




Low-angle electron diffraction



Differential phase contrast imaging









- STEM mode focused probe on sample, probe semi angle α .
- Beam deflected by Lorentz force.
- Segmented detector can then be used to map deflection by taking difference signals from opposite segments (quadrants or halves)










HF-etched FeRh planar sample as a function of temperature






Medipix pixelated detector can remove effects of diffraction contrast



HF-etched FeRh examined using the Medipix detector at 80°C






Cross-section of **55 nm FeRh** on **MgO** substrate


20°C \longrightarrow 200°C



200°C \longrightarrow 20°C




Planar TEM lamella of **FeRh** on **NiAl** on **GaAs** substrate at 150°C



1 μm



Biasing - 500 μ s pulses increasing from 100mV to 1.1V




↑
100 mV

1 μ m



Biasing - 500 μ s pulses increasing from 100mV to 1.1V




↑
400 mV

1 μ m




Biasing - 500 μ s pulses increasing from 100mV to 1.1V



1 μ m



Biasing - 500 μ s pulses increasing from 100mV to 1.1V




↑
800 mV

1 μ m



Biasing - 500 μ s pulses increasing from 100mV to 1.1V




↑
900 mV

1 μ m



Biasing - 500 μ s pulses increasing from 100mV to 1.1V




↑
1 V

1 μ m




Biasing - 500 μ s pulses increasing from 100mV to 1.1V





1 μ m

↑
1.1 V



- **Specialised TEM techniques** can be used to **visualise nanoscale magnetism** and provide fundamental insight for a range of applications
- Combining **electron holography** with *in situ* TEM and ETEM provided direct access to:
 - The effect of oxidation on vortex-state Fe_3O_4 grains
 - Thermomagnetic behaviour of vortex states
- Lorentz techniques like Fresnel, SAES, conventional and pixelated DPC imaging allow for:
 - Visualising the **dynamic nucleation and agglomeration** ferromagnetic domains during transition in FeRh
 - **DW and phase boundary motion** can be controlled by current pulsing and **Ir / Pd** doping of the FeRh films



Special Thanks to:

Adrian Muxworthy

**Imperial College
London**

**Wyn Williams
Lesleis Nagy**



**Stephen McVitie
Damien McGrouther
Sam McFadzean
Gary Paterson**



**University
of Glasgow**

Special Thanks to:

Takeshi Kasama
Thomas Hansen
Jakob Wagner
Jens Kling
Christian Damsgaard



Technical University of Denmark

Rafal Dunin-Borkowski
András Kovács
Chris Boothroyd
Vadim Migunov





Thank you